Mediloy M-Ti4

Ti100,0 [%]

CE 0197

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CE 0197

Instructions for use

Dental Ti-based metal-ceramic milling material, Type 4
Mediloy M-Ti4 is available as blanks.
Mediloy M-Ti4 complies with ISO 22674 and ISO 9693-1.

REF

Blanks, Ø 98 mm

Height (mm)	10	12	14	16
With shoulder	—	50571	50572	50573
Height (mm)	18	20	22	25
With shoulder	50574	50575	50576	50577

Alloy characteristics

According to ISO 22674 free of nickel, cadmium, beryllium and lead

Type (accord. to ISO 22674)	4
Melting point	°C 1668
Density	g/cm³ 4.5
Modulus of elasticity	GPa 120/115*
0.2% elongation limit ($R_{p,0.2}$) MPa	635/580*
Ductile yield (A_3)	% 20/18*
Hardness	(HV10) 225/235*
BEGO color code	8 (white)
Coefficient of thermal expansion (CTE)	
25 – 500 °C 10⁻⁶ K⁻¹	9.1 / 9.1*
(as machined / * simulated ceramic firings)	
Veneering ceramic	Ceramic with suitable CTE, e. g.: VITA Titankeramik
Oxidation firing	no oxide or control firing applicable!
Highest recommended firing temperature	840 °C
Heating rate	recommended max. 55 °C/min
Laser wire:	Titanwire Grade 2 Ø 0.35 mm (REF 50008)

Intended use: Mediloy M-Ti4 is indicated for the fabrication of dental restorations by milling.

Indication: Mediloy M-Ti4 is an unalloyed dental titanium for milling. It is suitable for the fabrication of crowns and bridges or metal-ceramic crowns and bridges, as well as implant prosthetics (e.g. abutments and bars).

Contraindications: No contraindications are known. However, unwanted biological reactions such as allergies to contents of the alloy or electrochemically based reactions may very rarely occur. In case of known incompatibilities and allergies to contents of the metallic material it should not be used.

Warnings: Metal dust is harmful to your health. When grinding and blasting use suitable air extraction system / ventilation at the workplace and breathing mask type FFP3-EN149! Titanium is highly reactive to oxygen when hot: fire and deflagration hazard! When using exhaust systems change filters regularly. Important! Fire hazard! Titanium is corroded by fluorides and hydrofluoric acid. Never use hydrofluoric acid for pickling, avoid contact to fluoride containing cleaning agents and/or dental care products.

Precautions: In case of occlusal or approximal contact with a different alloy electrochemically induced sensitivity may very rarely occur. Mediloy M-Ti4 has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of Mediloy M-Ti4 in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

Adverse reactions: No adverse reactions are known. Nevertheless, the rare case of occurrence of individual reactions against single components of Mediloy M-Ti4 cannot be excluded completely. In this case, the application of Mediloy M-Ti4 should not be continued.

Prescription device: Caution: US Federal law restricts this device to sale by or on the order of a licensed dentist.

Digital wax up: Minimum metal thickness (after grinding) 0.4 mm, implant prosthetics with screw channels 0.5 mm. Avoid sharp edges and corners. Framework should be anatomically reduced. Connectors should be modeled as strong and high as possible (height: min 3.5 mm, width: min. 2.5 mm).

Storage conditions: Dry.

Caution: Do not rework abutment interfaces!

Milling process: Mediloy M-Ti4 is compatible with following milling machines: Röders (e.g. RXD5, RXD5C). Please follow the instructions of the manufacturer of the devices. Clean thoroughly (steam clean or boil in aqua dest.).

Caution when removing blanks: The material may be hot or sharp-edged.

Use suitable hard metal cutters, e.g. with titanium coating. When using cooling/lubricating agents for milling, ensure complete final cleaning of surfaces!

Removal of restorations from blanks: Avoid dust formation! After cooling down, remove the restorations using e. g. a rotary cutter.

Finishing: Minimum metal thickness 0.4 mm. Use fine carbide tools suitable for Titanium. For grinding apply low pressure and work in one direction only. Clean thoroughly (steam clean or boil in aqua dest.).

Important: Use special tools for Titanium only. Keep these tools separate! A dark or discolored, slightly rough surface indicates incorrect finishing. Repeat finishing!

Passivation: Provide for sufficient passivation of a freshly prepared surface at different working stages by waiting for at least 10 min before resuming work.

Ceramics: For abutments and implant-supported, screw-retained bridges, ceramic firings must not be carried out! Use Titanium ceramics with suitable CTE according to ISO 9693-1. Always follow the ceramic manufacturer's instructions! Always blast the surface to be veneered with pure Korox® 110 (2 bar) thoroughly and with caution, allow passivation of the fresh surface for at least 10 min and then clean the framework thoroughly (steam clean or boil in aqua dest.). After cleaning, hold with artery clamps and refrain from touching.

Oxide firing: Do not process.

Ceramic firing: Please follow the in ceramics manufacturers' instructions! Remove ceramics only mechanically. Hydrofluoric acid (HF) strongly corrodes the metal frame.

Acrylic veneering: The corresponding instructions of the manufacturers must be followed when processing the veneering systems.

Final work: Carefully blast visible metal surfaces (e. g. inside crowns) with Korox® 50 at 2 bar and than blast-polish external surfaces with Perlblast®. This is of course dispensable for external surfaces when it has been completed before firing. Then continue polishing the external surfaces with low pressure and low speed (avoid overheating!) using polishers suitable for Titanium, after that finish-polish with BEGO-cobalt-chrome polishing paste (blue, REF 52310) and allow passivation of the fresh surface for at least 10 min. Clean thoroughly (steam clean or boil in aqua dest.).

Laser welding: If applicable use X-seam and filler material. Follow manufacturer's instructions for use and hazard notes of the laser welder devices. Clean thoroughly (steam clean or boil in aqua dest.).

No reuse of milled material: Already milled material (e. g. a bridge or bar) shall not be used for re-fabrication of dental restorations (e. g. by casting).

The blank may be used after several clampings until no more dental restorations can be milled from the remaining material.

Limit of Liability: Except where prohibited by law, BEGO Bremer Goldschlägerei Wilh. Herbst GmbH & Co. KG will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

Warranty: Whether given verbally, in writing or by practical instructions, our recommendations for use are based upon our own experience and trials and can be considered as standard values. Our products are subject to a constant further development. Therefore alterations in construction and composition are reserved.

US Labeling requirements: The device labeling meets the recommendations of FDA applicable guidance documents.

Any serious incident that has occurred in relation to Mediloy M-Ti4 should be reported to BEGO Bremer Goldschlägerei Wilh. Herbst GmbH & Co. KG and the competent authority.

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