Mediloy® Titan Grade 5 SLM

Ti89.8Al6.0V4.0Fe0 [%]

CE 0197

de Gebrauchsanweisung	2
en Instructions for use	3
bg Инструкции за употреба	4
cs Návod k použití	5
da Brugsanvisning	6
el Οδηγίες χρήσης	7
es Instrucciones de uso	8
et Kasutusjuhend	9
fi Käyttöohje	10
fr Notice d'utilisation	11
hr Upute za uporabu	12
hu Használati útmutató	13
it Istruzioni per l'uso	14
lt Naudojimo instrukcija	15
lv Lietošanas pamācība	16
nl Gebruiksaanwijzing	17
pl Instrukcja używania	18
pt Instruções de utilização	19
ro Instrucțiuni de utilizare	20
sk Návod na použitie	21
sl Navodila za uporabo	22
sv Bruksanvisning	23

Mediloy® Titan Grade 5 SLM

Ti89.8Al6.0V4.0Fe0 [%]

CE 0197

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CE 0197

Instructions for use

Titanium-based dental alloy for metal-ceramics, type 4

Grain size 20 – 63 µm

Mediloy® Titan Grade 5 SLM complies with ISO 22674 and ISO 9693.

REF 50525 - 2.5 kg

Alloy characteristics after heat treatment		900°C
According to ISO 22674 free of nickel, cadmium, beryllium and lead		
Type (according to ISO 22674)		4
Density	g/cm ³	4.3
Solidus, liquidus temperature	°C	1604, 1660
Modulus of elasticity	GPa	120/115*
0.2% elongation limit (R _{p0.2})	MPa	875/880*
Ductile yield (A ₅)	%	11/10*
Vickers hardness	HV10	320/320*
BEGO colour code		8 (white)
Coefficient of thermal expansion (CTE)		
25 – 500°C	10 ⁻⁶ K ⁻¹	10.3/10.0*
Highest recommended firing temperature	°C	840
Heating rate		recommended max. 55°C/min
Veneering ceramic		Ceramics with compatible CTE value, e.g.: VITA LUMEX AC
Oxide firing		Not recommended
Laser wire		Titandraht Grade 2 (REF 50008)

*Measured value after simulated ceramic firings

Intended purpose: The device is suitable for producing dental restorations, such as crowns, bridges and metal-ceramic frames with the SLM method (selective laser melting).

Indication: The device is intended for treating missing hard tissue (teeth).

Contraindication: No contraindications are known.

Clinical benefit: Artificial replacement for lost hard tissue (teeth), for restoring masticatory function (aesthetic and functional).

Warnings: Metal dust is harmful to health. Avoid dust formation! The opening of packages, filling of powders, grinding and blasting of dental restorations should be performed carefully and using an appropriate extraction system. Respiratory protection of type FFP3-EN149, protective goggles with side protection (DIN EN 166), protective gloves (made of butyl rubber or nitrile rubber, category III, EN 374) and ESD-certified safety shoes are recommended.

When hot, titanium powder reacts to oxygen by annealing. Risk of burns and risk of deflagration! When using extraction systems, replace the paper filter regularly. Caution! Otherwise risk of fire! Titanium is corroded by fluorides and hydrofluoric acids: never etch with hydrofluoric acid, do not allow to come into contact with cleaning agents/dental hygiene products with fluoride contents.

In the event of contact with eyes, rinse with plenty of water. In the event of skin contact, wash with water and soap. If irritation persists, seek a physician's care.

Collect any spilled amounts mechanically with a damp cloth (isopropanol) or with an ATEX vacuum cleaner, and dispose of in accordance with local and national statutory regulations. Metal powders are combustible. Remove all sources of ignition. Suitable extinguishing media: special powders against metal fires, sand. Pay attention to safety data sheet!

Precautions: In the case of approximal or occlusal contact with other metals, electrochemically related numbness may occur in very rare cases. The device has not been assessed for safety and compatibility in the MR environment. It has not been tested for temperature increase, migration or image artefacts in the MR environment. Safety in the MR environment is unknown. Examination of a patient with this device can lead to injuries to the patient. Faulty restorations must not be inserted into the mouth.

Patient group: Objects made of Mediloy® Titan Grade 5 SLM can be used regardless of the age and gender of the patient. The alloy should not be used in the event of known intolerances or known allergies to alloy components.

Side effects: No side effects are known. However, the possibility of individual reactions to components in isolated cases cannot be ruled out. In such cases, the device should not be used.

Digital wax-up: Wax-up is performed using suitable CAD software under consideration of dental technology regulations. Wall strengths after finishing of at least 0.4 mm, avoid sharp corners and edges. Veneer frames to be designed in anatomically reduced form. Allow connectors to be as strong and high as possible (height: min. 3.5 mm, width: min. 2.5 mm).

Work steps in the manufacturing centre: for equipment-specific work steps and settings, the device manufacturer's specifications must be adhered to! Please follow the instructions for use and safety instructions of the equipment manufacturers!

Storage conditions: Keep dry in tightly closed container.

SLM method: Prevent the formation of dust when opening the packaging, during transport and when filling the powder into the SLM system. Use an SLM system with suitable laser (e.g. Ytterbium fibre laser or Nd:YAG laser (wavelength approximately 1060 – 1100 nm)) with the following settings: **Laser output 190 W**, laser beam diameter 0.08 mm (e.g. Renishaw), track spacing 0.07 mm, powder coating thickness 0.03 mm, scanning speed 1,600 mm/s, energy density 1.5 Joule/mm², platform heating 100°C. The values are exemplary for common devices. Therefore, the current operating instructions provided by the device manufacturer must be observed for the individual settings for titanium powder alloying. A suitable inert gas for titanium powder alloying must be used.

Cleaning components: Remove the production platform from the SLM system, pour out unused powder and vacuum-clean the platform (objects) using an ATEX vacuum cleaner to remove all residue powder.

Reuse: If unmelted powder is to be reused, it should be sifted beforehand using an ATEX ultrasound sieve (75 µm) and larger grain sizes are ejected. Also make sure that unmelted powder for reuse does not become contaminated by foreign particles.

Heat treatment: The removable part of the production platform with the produced objects is placed in a suitable kiln (e.g. hot wall retort kiln, argon or vacuum) at room temperature. The temperature is increased to 900°C within 30 minutes, and then held for 60 minutes. The temperature is then decreased to room temperature within 150 minutes. The platform can then be removed from the kiln for further processing.

Separating the restorations from the plate: Prevent the formation of dust! After the heat treatment and cooling of the platform, remove the restorations using a band saw, rotary instruments or forceps, for example. Also remove the remaining supports using forceps.

No re-use of laser-sintered material: Material (e.g. a bridge or residue of melted powders) which has already been melted by SLM, must not be reused for the manufacture of a new restoration (e.g. by casting). The alloy characteristics can only be guaranteed for new material.

Finishing: Use cross-cut carbide burs which are suitable for titanium, and work with low contact pressure in one direction.

Caution: Only use special burs for titanium machining! Keep separate! A dark or discoloured, slightly rough surface indicates insufficient machining. Repeat finishing.

Passivation: Sufficient passivation must always be guaranteed after various working steps. To do so, allow the freshly machined surfaces to sit for at least 10 minutes.

Ceramic veneering: Use veneering ceramics with a suitable CTE (ISO 9693), observe the instructions for use of the respective ceramics manufacturer. Before ceramic veneering, it is imperative that the surfaced being veneered first be carefully and thoroughly blasted (clean Korox® 110 at 2 bar), allow to sit for at least 10 minutes for passivation. Thorough cleaning with a steam blaster, ultrasound bath or by boiling in distilled water is required. After this step, do not touch the surfaces again with your hands. Use arterial clamps or similar. Ensure that the frames are supported appropriately during firing.

Oxide firing: Do not carry out.

Firing: Please follow the instructions of the ceramics manufacturer. Only remove ceramics mechanically. Hydrofluoric acid (HF) heavily corrodes the metal frame.

Composite veneers: The respective manufacturers' instructions must be heeded when working with composite veneering materials.

Concluding work: Carefully blast visible metal surfaces (e.g. inner surfaces of crowns) with Korox® 50 at 2 bar, and blast polish outer surfaces with Perlblast®. Then machine with polishers which are suitable for titanium, with low contact pressure and at a low speed (avoid overheating!), then polish with BEGO blue polishing compound (REF 52310) and allow to sit for at least 10 minutes for passivation. Thorough cleaning of the surface with a steam blaster or by boiling in distilled water.

Laser welding: If possible, work with X sutures and filler alloy (e.g. Titandraht Grade 2; REF 50008). Please follow the instructions for use and safety instructions of the equipment manufacturer!

Cleaning/disinfection: The restorations are manufactured non-sterile and must be dismantled into individual parts, cleaned and disinfected before insertion into the patient's mouth.

Note: Each prosthetics component must only be used for its designated purpose. After machining, the prosthetic work must be cleaned with a steam blaster, in an ultrasound bath or by boiling in distilled water, as appropriate, and disinfected before insertion into the patient, as applicable.

Disinfection: When selecting the disinfectant, you must make sure that it is suitable for cleaning and disinfection of dental restorations, and is compatible with the materials of the restorations being cleaned and disinfected. The instructions of the disinfectant manufacturer must be observed.

Warranty: Application-related recommendations provided by us, whether given verbally, in writing or by way of practical instructions, are based on our own experience and tests and may therefore only be regarded as general guidelines. Our products are subject to continuous development. We thus reserve the right to make modifications in design, appearance and materials without notice. Please report any serious incident that has occurred in relation to the device to BEGO Bremer Goldschlägerei Wilh. Herbst GmbH & Co. KG and the competent authority.

Information on disposal: Procedure for waste treatment

Device: The assignment of a waste key number as per the European Waste Catalogue Ordinance (AVV) must be carried out in consultation with the regional waste disposal contractor. Do not dispose of with household waste.

Packaging: Packaging must be completely emptied and disposed of properly in accordance with the statutory regulations. Packaging that is not fully emptied must be disposed of in coordination with the regional waste disposal contractor.

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